

Used Vision Tester, Refurbished
 +50 shipping
\$499

Vision Tester Features:

- + It is very simple & easy to operate
- + It has pointer hole on side, to point to slides
- + It is durable and sturdy
- + It is easy to service and repair
- + It has near & far vision switch
- + Slides can easily be replaced
- + It is compact & smaller than most vision testers
- + It is easy to train nurses with this machine
- + It has all the standard features that most nurses want
- + It is priced under \$1,000, whereas new machines cost \$2-4,000











Titmus 2S - Vision Tester
 The most reliable vision tester
 that Titmus has made in the past 40 years.

Professional vision testing slides
 Titmus models 2S, V2 V4

1. (Test No. 1, SA-APS-1) Allen Picture Test measures binocularity and acuity. Line 1 - Binocularity - tests performance of both eyes together. Lines 2,3,4 - Acuity - right and left eyes, calibrated at 20/30, 20/40 and 20/50.
2. (Test No. 2, SAS-2) meaasures Acuity in right and left eyes with Conventional Tumbling "E" symbol. Patient states which way the "E" points in each square.
3. (Test No. 3, SPS-3) measures Heterophorias. Muscle imbalance - vertical and lateral - is tested at predetermind pass/fail levels. This test is administered first at far, then at near.
4. (Test No. 4, SCD-1) is used to test color perception, and consists of six accurately and authentically reproduced Ishihara Pseudo-Isochromatic plates. Detects color deficiency, but does not classify.
5. (Test No. 5, SA-BRL-4) measure visual acuity. Results are calibrated at 20/20, 20/30, 20/40, 20/50, 20/70, 20/100 and 20/200. Tests can be administerd binocularly or with either eye occluded.
6. (Test No. 6, SLP-1) measures lateral phoria. Desirable phoria scores are associated with quick, accurate eye fixations and precise, effortless seeing.
7. (Test No. 7, SVP-1) measures vertical balance between the two eyes. Good vertical alignment is an important factor in comfortable, efficient, two-eyed vision and is related to efficiency in most occupations.
8. (Test No. 8, SA-BRL-2A) measures visual acuity. Visual acuity test results are calibrated at 20/20, 20/30, 20/40, 20/50, 20/70 and 20/100. Tests can be administerd binocularly or with either eye occluded.

Slides for Professional Model

 Test No. 1, SA-APS-1 Allen Picture Test Line 1 - Binocularity - Tests performance of both eyes together. Lines 2, 3, 4 - Acuity - Right and Left Eye Calibrated at 20/30, 20/40 and 20/50.	 Test No. 5, SA-BRL-4 - Visual Acuity Results calibrated at 20/20, 20/30, 20/40, 20/50, 20/70, 20/100, and 20/200. Tests can be administered binocularly or with either eye occluded.
 Test No. 2, SAS-2 - Acuity Right and Left Eye Conventional Tumbling "E" symbol. Child states which way the "E" points in each square. Acuity levels of 20/20, 20/30, 20/40 and 20/50.	 Test No. 6, SLP-1 - Lateral Phoria Measures Esophoria or Exophoria in 1 prism diopter steps.
 Test No. 3, SPS-3 - Heterophorias Muscle imbalance - vertical and lateral - is tested at predetermined pass/fail levels. This test is administered first at far, then at near.	 Test No. 7, SVP-1 - Vertical Phoria Measures right or left Hyperphorias in one-half prism diopter steps.
 Test No. 4, SCD-1 - Color Perception Consists of six accurately and authentically reproduced Ishihara Pseudo-Isochromatic plates. Detects color deficiency but does not classify.	 Test No. 8, SA-BRL-2A - Visual Acuity Results calibrated at 20/20, 20/30, 20/40, 20/50, 20/70, and 20/100. Tests can be administered binocularly or with either eye occluded.

Call us to confirm price and order
631-423-0797
 Services@HopkinsCalibration.com

Slides for Professional Model



Test No. 1, SA-APS-1 Allen Picture Test

Line 1 – Binocularity – Tests performance of both eyes together. Lines 2, 3, 4 – Acuity – Right and Left Eye
 Calibrated at 20/30, 20/40 and 20/50.



Test No. 5, SA-BRL-4 – Visual Acuity

Results calibrated at 20/20, 20/30, 20/40, 20/50, 20/70, 20/100, and 20/200.
 Tests can be administered binocularly or with either eye occluded.



Test No. 2, SAS-2 – Acuity Right and Left Eye

Conventional Tumbling “E” symbol.
 Child states which way the “E” points in each square.
 Acuity levels of 20/20, 20/30, 20/40 and 20/50.



Test No. 6, SLP-1 – Lateral Phoria

Measures Esophoria or Exophoria in 1 prism diopter steps.



Test No. 3, SPS-3 – Heterophorias

Muscle imbalance – vertical and lateral – is tested at predetermined pass/fail levels.
 This test is administered first at far, then at near.



Test No. 7, SVP-1 – Vertical Phoria

Measures right or left Hyperphorias in one-half prism diopter steps.



Test No. 4, SCD-1 – Color Perception

Consists of six accurately and authentically reproduced Ishihara Pseudo-Isochromatic plates.
 Detects color deficiency but does not classify.



Test No. 8, SA-BRL-2A – Visual Acuity

Results calibrated at 20/20, 20/30, 20/40, 20/50, 20/70, and 20/100.
 Tests can be administered binocularly or with either eye occluded.